

Oxford Cryosystems serial line communication protocols[home](#) | [700 and 800 series](#) | [GM cryocoolers](#)700 and 800 series: [home](#) | [protocols](#)[Cryostream](#) | [Cobra](#) | [N-HeliX](#) | [PheniX](#) | [PheniX-FL](#) | [Chimera](#) | [Smartstream](#)[status packets](#) | [command packets](#)

Cryostream command packets

All 700 and 800 series controllers respond to the same set of commands, with minor variations noted below. The structure of a valid command packet is as follows:

```
char    Size;          /* The variable size of the command packet */
char    Id;            /* The command packet identifier */
char    *Params;       /* Possible variable number of parameters */
```

The Size and Id parameters will match one of those listed in the table below. If a command is unrecognised (Id invalid or Size inappropriate), illegal (parameter out of range) or inappropriate (eg the machine has shutdown), then it is ignored. No command acknowledgement is issued, instead consult the next status packet for evidence that the command has been received.

Note: this information is available in machine-readable form in [Cryostream.xml](#).

List of commands

Command	Size	Id	Params	Notes
Restart	2	10	none	
Ramp	6	11	HIBYTE (RampRate), LOBYTE (RampRate), HIBYTE (TargetTemp), LOBYTE (TargetTemp)	
Plat	4	12	HIBYTE (Duration), LOBYTE (Duration)	
Hold	2	13	none	
Cool	4	14	HIBYTE (TargetTemp), LOBYTE (TargetTemp)	
End	2	15	none	
Purge	2	16	none	
Pause	2	17	none	
Resume	2	18	none	
Stop	2	19	none	
Turbo	3	20	TurboOn	
SetFormat	3	40	Format	Switch to extended status packets

Parameter units and ranges

Parameter	Units	Min	Max
RampRate	K/hour	1	360
TargetTemp	cK	[Cryostream] 8000	[Cryostream] 40000 or 50000 (Plus)
Duration	minutes	1	1440 (=24 hours)
TurboOn	none	0	1
Format	none	0 (old format)	1 (extended format)

Example commands

In most of the above cases no parameters are required, and thus the command packet is a simple two-byte

one in which Size=2. For example, a command packet to stop the cooler would be created as follows:

```
/*
 * The Stop command packet, size = 2, id = 19
 */
char buf[2] = { 2, 19 };    /* Create a Stop command packet */
```

For those commands requiring parameters, the Params[] array takes various forms, illustrated by the following examples.

The Turbo command takes a single 1-byte parameter:

```
/*
 * The Turbo command packet, size = 3, id = 20
 * The Params[] array consists of a single char taking the value either
 * 0 (switch Turbo off) or 1 (switch Turbo on)
 */
char buf[3] = { 3, 20, 1 };    /* Switches Turbo on */
```

The Ramp, Plat and Cool commands are a little more complicated, because the parameters which are passed are shorts - 16 bit integers. These need to be assembled as illustrated in the following examples, which use the standard [HIBYTE](#) and [LOBYTE](#) macros reproduced below.

```
/*
 * The Plat command packet, size = 4, id = 12
 * The Params[] array consists of a short containing the duration of
 * the
 * Plat in minutues
 */
char buf[4] = { 4, 12, HIBYTE(720), LOBYTE(720) }; /* 720 minute
plateau */

/*
 * The Cool command packet, size = 4, id = 14
 * The Params[] array consists of a short containing the end
 * temperature in centi-Kelvin
 */
char buf[4] = { 4, 14, HIBYTE(17000), LOBYTE(17000) }; /* Cool to 170 K
*/

/*
 * The Ramp command packet, size = 6, id = 11
 * The Params[] array consists of a short containing desired ramp rate
 * in K/hour,
 * followed by a short containing the end temperature in centi-Kelvin
 */
char buf[6] = { 6, 11, HIBYTE(120), LOBYTE(120), /* Rate = 120 K/hour
*/
                HIBYTE(25050), LOBYTE(25050), /* Final temp = 250.5
K */ };
```

HIBYTE

```
#define HIBYTE(w) (((unsigned char)((((unsigned short)(w) >> 8) & 0xFF))
```

LOBYTE

```
#define LOBYTE(w) ((unsigned char)(w))
```